
Product Datasheet

CF150 Antibody (Center) (orb1788102)

Description

Affinity Purified Rabbit Polyclonal Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	266-295 aa
Target	This CF150 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 266-295 amino acids from the Central region of human CF150.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	58814 Da
Uniprot ID	Q8N884
NCBI	NP_612450.2
Dilution Range	WB: 1:1000, FC: 1:25
Expiration Date	12 months from date of receipt.