
Product Datasheet

ENOB Antibody (Center) (orb1788251)

Description

Purified Rabbit Polyclonal Antibody (Pab)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, IHC-P, WB

Immunogen

237-264 aa

Target

This ENOB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 237-264 amino acids from the Central region of human ENOB.

Form/Appearance

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles

Note

For research use only

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

46987 Da

Uniprot ID[P13929](#)**NCBI**[NP_001967.3](#), [NP_001180432.1](#),
[NP_443739.3](#)**Dilution Range**

WB: 1:1000, IHC-P: 1:50~100, FC: 1:10~50

Expiration Date

12 months from date of receipt.