

Product Datasheet

Functional ADAM17 (human) Antibody, mAb (recombinant) (blocking) (preservative free) (orb1788485)

Description	Functional ADAM17 (human) Antibody, mAb (recombinant) (blocking) (preservative free)
Reactivity	Human
Conjugation	Unconjugated
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P78536
Expiration Date	12 months from date of receipt.