

Product Datasheet

Functional PEDF (human) Antibody, mAb (recombinant) (orb1788487)

Description	Functional PEDF (human) Antibody, mAb (recombinant)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P36955
Expiration Date	12 months from date of receipt.