

Product Datasheet

MUC5AC Antibody (orb1789793)

Description	Mouse Monoclonal Antibody [Clone MUC5AC/917 + 45M1]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P98088
Expiration Date	12 months from date of receipt.