

Product Datasheet

HLA-B Antibody (orb1789953)

Description	Mouse Monoclonal Antibody [Clone CATA-1]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Uniprot ID	P01889
Expiration Date	12 months from date of receipt.