

Product Datasheet

TAG-72 / CA72.4 (Tumor-Associated Glycoprotein) Antibody - With BSA and Azide (orb1790034)

Description	Mouse Monoclonal Antibody [Clone CA72/145]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IP, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.