

Product Datasheet

Bromodeoxyuridine (BrdU) (Proliferation Marker) Antibody - With BSA and Azide (orb1790089)

Description	Mouse Monoclonal Antibody [Clone BRD469 + BRD494 + BRD.3]
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles
Note	For research use only
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.