

Product Datasheet

Human Alkaline Phosphatase protein (orb180168)

Description	Recombinant Human Alkaline Phosphatase is produced by our Mammalian expression system and the target...
Endotoxins	Less than 0.1 EU/ug
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM Tris-Hcl, 150mM NaCl, 0.2mM CaCl ₂ , 10% Glycerol, pH 7.5
Storage	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Tag	C-6His
Note	For research use only
Uniprot ID	P05187
Expiration Date	6 months from date of receipt.