

## Product Datasheet

### Human ABCB5 protein (orb180175)

<b>Description</b>	Recombinant Human ATP-binding cassette sub-family B member 5 is produced by our E.coli expression sy...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	N-Trx
<b>Note</b>	For research use only
<b>Application notes</b>	Trx tag at the N-terminus.
<b>Uniprot ID</b>	<a href="#">Q2M3G0</a>
<b>Expiration Date</b>	6 months from date of receipt.