

## Product Datasheet

### Human CAMK1D protein (orb180197)

<b>Description</b>	Recombinant Human Calcium/Calmodulin-Dependent Protein Kinase Type ID is produced by our E.coli expr...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Tag</b>	N-GST
<b>Note</b>	For research use only
<b>Application notes</b>	GST tag at the N-terminus.
<b>Uniprot ID</b>	<a href="#">Q8IU85</a>
<b>Expiration Date</b>	6 months from date of receipt.