

Product Datasheet

Human CLIC3 protein (orb180228)

Description	Recombinant Human Chloride intracellular channel protein 3 is produced by our E.coli expression syst...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 10mM Tris-Hcl, 0.1% Triton X-100, pH 8.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	O95833
Expiration Date	6 months from date of receipt.