
Product Datasheet

C Reactive Protein/CRP Antibody (orb18025)

Description

C Reactive Protein/CRP Antibody

Species/Host

Rabbit

Reactivity


Human

Conjugation

Unconjugated

Tested

WB

Applications


Immunohistochemical analysis of paraffin...

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human C Reactive Protein (209-224aa KYEVQGEVFTKPLWP).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.


Isotype

Rabbit IgG

Clonality

Polyclonal

MW

25039 MW

Uniprot ID
[P02741](#)
Sensitivity

10pg/ml

Expiration Date

12 months from date of receipt.