

Product Datasheet

Human GABARAP (N-His tag) protein (orb180276)

Description	Recombinant Human GABA(A) receptor-Associated Protein is produced by our E.coli expression system an...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB, 150mM NaCl, 20% Glycerol, pH 7.0.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His&C-Fc
Note	For research use only
Application notes	6His tag at the N-terminus, Fc tag at the C-terminus.
Uniprot ID	Q6IAW1
Expiration Date	6 months from date of receipt.