

Product Datasheet

Human UBE2M protein (orb180352)

Description	Recombinant Human Ubiquitin-conjugating Enzyme E2M is produced by our E.coli expression system and t...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 50mM HEPES, pH 7.5, 2mM DTT, 150mM NaCl, 10% Glycerol.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	tag free
Note	For research use only
Uniprot ID	P61081
Expiration Date	6 months from date of receipt.