

## Product Datasheet

### Human SORD protein (orb180413)

<b>Description</b>	Recombinant Human Sorbitol Dehydrogenase is produced by our Mammalian expression system and the targ...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 200mM NaCl, 5mM DTT, 20% Glycerol, pH 8.0.
<b>Storage</b>	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Application notes</b>	6His tag at the C-terminus.
<b>Uniprot ID</b>	AAH21085.1
<b>Expiration Date</b>	6 months from date of receipt.