

Product Datasheet

Human UBE2B protein (orb180450)

Description	Recombinant Human Ubiquitin-Conjugating Enzyme E2 B is produced by our Mammalian expression system a...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM PB,150mM NaCl, pH 7.4.
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	C-His
Note	For research use only
Application notes	6His tag at the C-terminus.
Uniprot ID	P63146
Expiration Date	6 months from date of receipt.