

Product Datasheet

Human UBA5 protein (orb180454)

Description	Recombinant Human Ubiquitin-fold Modifier 1 Activating Enzyme is produced by our E.coli expression s...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Supplied as a 0.2 um filtered solution of 20mM TrisHCl, 50mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0 .
Storage	Store at -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Tag	N-His
Note	For research use only
Application notes	6His tag at the N-terminus.
Uniprot ID	Q9GZZ9
Expiration Date	6 months from date of receipt.