



Product Datasheet

Anti-GAPDH Antibody (HRP) (orb1804652)

Description	0
Species/Host	Rabbit
Reactivity	Gallus, Human, Monkey, Mouse, Rat, Zebrafish
Conjugation	HRP
Tested Applications	WB
Immunogen	HRP-GAPDH antibody was raised against a synthetic peptide containing 16 amino acids near the carboxy terminus of GAPDH.
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	HRP-GAPDH antibody can be stored at 4°C for three months and -20°C, stable for up to one year.
Note	For research use only
Application notes	Western blot, 0.02-0.1 µg/ml
Isotype	IgG
Clonality	Polyclonal
MW	36053 Da
Uniprot ID	P04406
Expiration Date	12 months from date of receipt.



WB analysis of GAPDH using anti-GAPDH an...



WB analysis of GAPDH using anti-GAPDH an...