

Product Datasheet

Retrorsine Antibody (orb1804724)

Description	Retrorsine antibody
Species/Host	Rabbit
Reactivity	Animal
Conjugation	Unconjugated
Tested Applications	ELISA, RIA
Immunogen	Retrorsin-hemi-CDI-BSA or BTG
Target	Retrorsine
Preservatives	10mM sodium azide. PBS, pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	1kDA
Expiration Date	12 months from date of receipt.