

Product Datasheet

Nanosilver Antibody (orb1804743)

Description	nanosilver antibody
Species/Host	Rabbit
Reactivity	Animal
Conjugation	Unconjugated
Tested Applications	ELISA, RIA
Immunogen	Nanosilver-3MBA-BSA or KLH
Target	nanosilver
Preservatives	10mM sodium azide. PBS, pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	1kDA
Expiration Date	12 months from date of receipt.