

Product Datasheet

T-2 toxin Antibody (orb1804777)

| | |
|----------------------------|---------------------------------|
| Description | T-2 toxin antibody |
| Species/Host | Rabbit |
| Reactivity | Animal |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, RIA |
| Immunogen | T2-HS-KLH |
| Target | T-2 toxin |
| Preservatives | 10mM sodium azide. PBS, pH 7.4 |
| Form/Appearance | Liquid |
| Storage | -20°C |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 1kDA |
| Expiration Date | 12 months from date of receipt. |