

Product Datasheet

Fumonisin Antibody (orb1804781)

| | |
|----------------------------|---------------------------------|
| Description | Fumonisin antibody |
| Species/Host | Rabbit |
| Reactivity | Animal |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, RIA |
| Immunogen | FB1-EDC-LPH |
| Target | Fumonisin |
| Preservatives | 10mM sodium azide. PBS, pH 7.4 |
| Form/Appearance | Liquid |
| Storage | -20°C |
| Note | For research use only |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 1kDA |
| Expiration Date | 12 months from date of receipt. |