

## Product Datasheet

### Nitrofurans Antibody (orb1804785)

<b>Description</b>	Nitrofurans antibodies: 1.NPAOZ 2.SEM 3.AMOZ 4.AOZ 5.AHD
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Animal
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, RIA
<b>Target</b>	Nitrofurans : 1.NPAOZ 2.SEM 3.AMOZ 4.AOZ 5.AHD
<b>Preservatives</b>	10mM sodium azide. PBS, pH 7.4
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	1kDA
<b>Expiration Date</b>	12 months from date of receipt.