

Product Datasheet

Spectinomycin Antibody (orb1804799)

Description	Spectinomycin antibody
Species/Host	Rabbit
Reactivity	Animal
Conjugation	Unconjugated
Tested Applications	ELISA, RIA
Immunogen	Spectinomycin-CMO-BTG
Target	Spectinomycin
Preservatives	10mM sodium azide. PBS, pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	1kDA
Expiration Date	12 months from date of receipt.