

## Product Datasheet

### Streptomycine Antibody (orb1804804)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | Streptomycine antibody          |
| <b>Species/Host</b>        | Rabbit                          |
| <b>Reactivity</b>          | Animal                          |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | ELISA, RIA                      |
| <b>Immunogen</b>           | Streptomycin-KLH                |
| <b>Target</b>              | Streptomycine                   |
| <b>Preservatives</b>       | 10mM sodium azide. PBS, pH 7.4  |
| <b>Form/Appearance</b>     | Liquid                          |
| <b>Storage</b>             | -20°C                           |
| <b>Note</b>                | For research use only           |
| <b>Isotype</b>             | Rabbit IgG                      |
| <b>Clonality</b>           | Polyclonal                      |
| <b>MW</b>                  | 1kDA                            |
| <b>Expiration Date</b>     | 12 months from date of receipt. |