

Product Datasheet

Diethylstilbestrol Antibody (orb1804832)

Description	Diethylstilbestrol antibody
Species/Host	Rabbit
Reactivity	Animal
Conjugation	Unconjugated
Tested Applications	ELISA, RIA
Immunogen	diethylstilbestrol-CPE-BSA
Target	Diethylstilbestrol
Preservatives	10mM sodium azide. PBS, pH 7.4
Form/Appearance	Liquid
Storage	-20°C
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	1kDA
Expiration Date	12 months from date of receipt.