

## Product Datasheet

# SARS-CoV-2 (Omicron BA.2) Spike S1, His Tag (orb1805231)

<b>Description</b>	SARS-CoV-2 (Omicron BA.2) Spike S1, His Tag is expressed from HEK293. It contains AA Ser 13-Arg 682 ...
<b>Species/Host</b>	Human
<b>Reactivity</b>	Virus
<b>Endotoxins</b>	< 1 EU/mg
<b>Tested Applications</b>	ELISA, Kinetics
<b>Target</b>	Spike S1
<b>Preservatives</b>	1×PBS, pH 6.0
<b>Storage</b>	-80°C
<b>Tag</b>	His
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Ser 13 - Arg 682
<b>Isotype</b>	His
<b>Purity</b>	> 85%
<b>MW</b>	76.75 kDa
<b>Uniprot ID</b>	<a href="#">PODTC2</a>
<b>Expiration Date</b>	6 months from date of receipt.