

Product Datasheet

SARS-CoV-2 (Omicron B.1.1.529) Spike S1, Avi-His Tag (orb1805235)

Description	SARS-CoV-2 (Omicron B.1.1.529) Spike S1, Avi-His Tag is expressed from HEK293. It contains AA Ser 13...
Species/Host	Human
Reactivity	Virus
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Spike S1
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Ser 13 - Arg 682
Isotype	His
Purity	> 85%
MW	77.98 kDa
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.