

Product Datasheet

SARS-CoV-2 (Omicron B.1.1.529) Spike S1 (Q493K), Fc Tag (orb1805236)

Description SARS-CoV-2 (Omicron B.1.1.529) Spike S1 (Q493K), Fc Tag is expressed from HEK293. It contains AA Ser...

Species/Host Human

Reactivity Virus

Endotoxins < 1 EU/mg

Tested Applications ELISA, Kinetics

Target Spike S1

Preservatives 1×PBS, pH 6.0

Storage -80°C

Tag Fc

Note For research use only

Protein Sequence Ser 13 - Arg 682

Isotype Fc

Purity > 85%

MW 202.78 kDa

Uniprot ID [PODTC2](#)

Expiration Date 6 months from date of receipt.