

Product Datasheet

MonkeypoxVirus (strain Zaire-96-I-16) A30L, His Tag (orb1805271)

| | |
|----------------------------|---|
| Description | MonkeypoxVirus (strain Zaire-96-I-16) A30L, His Tag is expressed from HEK293. It contains AA Ile 26-... |
| Species/Host | Human |
| Reactivity | Virus |
| Endotoxins | < 1 EU/mg |
| Tested Applications | ELISA, Kinetics |
| Target | A30L |
| Preservatives | 1×PBS, pH 6.0 |
| Storage | -80°C |
| Tag | His |
| Note | For research use only |
| Protein Sequence | Ile 26-Leu 146 |
| Isotype | His |
| Purity | > 85% |
| MW | 15.18 kDa |
| Uniprot ID | Q8V4U9 |
| Expiration Date | 6 months from date of receipt. |