

Product Datasheet

Human SFRP1, His Tag (orb1805307)

Description	Human SFRP1, His Tag is expressed from HEK293. It contains AA Met 1-Lys 314 (Accession No. Q8N474).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	SFRP1
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Met 1-Lys 314
Isotype	His
Purity	> 85%
MW	34 kDa
Uniprot ID	Q8N474
Expiration Date	6 months from date of receipt.