

Product Datasheet

Human DKK3, His Tag (orb1805315)

Description	Human DKK3, His Tag is expressed from HEK293. It contains AA Ala22-Ile350 (Accession No. Q9UBP4).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	DKK3
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Ala22-Ile350
Isotype	His
Purity	> 85%
MW	37.22 kDa
Uniprot ID	Q9UBP4
Expiration Date	6 months from date of receipt.