

Product Datasheet

Human RNF43, Fc Tag (orb1805319)

Description	Human RNF43, Fc Tag is expressed from HEK293. It contains AA Met1-Tyr197 (Accession No. Q68DV7).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	RNF43
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Met1-Tyr197
Isotype	Fc
Purity	> 85%
MW	45.8 kDa
Uniprot ID	Q68DV7
Expiration Date	6 months from date of receipt.