

Product Datasheet

Human TNFRSF12A / TWEAKR / CD266, Fc Tag (orb1805320)

| | |
|----------------------------|--|
| Description | Human TNFRSF12A / TWEAKR / CD266, Fc Tag is expressed from HEK293. It contains AA 53 (Accession No. ...) |
| Species/Host | Human |
| Reactivity | Human |
| Endotoxins | < 1 EU/mg |
| Tested Applications | ELISA, Kinetics |
| Target | TNFRSF12A / TWEAKR / CD266 |
| Preservatives | 1×PBS, pH 6.0 |
| Storage | -80°C |
| Tag | Fc |
| Note | For research use only |
| Protein Sequence | Glu 28 - Pro 80 |
| Isotype | Fc |
| Purity | > 85% |
| MW | 63.36 kDa |
| Uniprot ID | Q9NP84 |
| Expiration Date | 6 months from date of receipt. |