

Product Datasheet

Human TPBG, His Tag (orb1805323)

Description	Human TPBG, His Tag is expressed from HEK293. It contains AA 324 (Accession No. Q13641).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	TPBG
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Ser 32 - Ser 355
Isotype	His
Purity	> 85%
MW	36.42 kDa
Uniprot ID	Q13641
Expiration Date	6 months from date of receipt.