

Product Datasheet

Human TLR4 / CD284, Fc Tag (orb1805324)

Description	Human TLR4 / CD284, Fc Tag is expressed from HEK293. It contains AA 608 (Accession No. O00206-1).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	TLR4 / CD284
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Glu 24 - Lys 631
Isotype	Fc
Purity	> 85%
MW	190.5 kDa
Uniprot ID	O00206
Expiration Date	6 months from date of receipt.