

Product Datasheet

Human Transferrin receptor / CD71, Fc Tag (orb1805326)

Description	Human Transferrin receptor / CD71, Fc Tag is expressed from HEK293. It contains AA 697 (Accession No...
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Transferrin receptor / CD71
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Arg 105 - Phe 801
Isotype	Fc
Purity	> 85%
MW	206.34 kDa
Uniprot ID	Q9UP52
Expiration Date	6 months from date of receipt.