

Product Datasheet

Human PTCRA, His Tag (orb1805333)

Description	Human PTCRA, His Tag is expressed from HEK293. It contains AA 123 (Accession No. Q6ISU1-1).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	PTCRA
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Thr 24 - Gly 146
Isotype	His
Purity	> 85%
MW	13.86 kDa
Uniprot ID	Q6ISU1
Expiration Date	6 months from date of receipt.