

Product Datasheet

Human PRLR / Prolactin Receptor, His Tag (orb1805334)

Description	Human PRLR / Prolactin Receptor, His Tag is expressed from HEK293. It contains AA 210 (Accession No....
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	PRLR / Prolactin Receptor
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Gln 25 - Asp 234
Isotype	His
Purity	> 85%
MW	25.51 kDa
Uniprot ID	P16471
Expiration Date	6 months from date of receipt.