

Product Datasheet

Human TNFRSF4 / OX40 / CD134, Fc Tag (orb1805338)

Description	Human TNFRSF4 / OX40 / CD134, Fc Tag is expressed from HEK293. It contains AA 186 (Accession No. P43...
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	TNFRSF4 / OX40 / CD134
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Leu 29 - Ala 214
Isotype	Fc
Purity	> 85%
MW	92.12 kDa
Uniprot ID	P43489
Expiration Date	6 months from date of receipt.