

Product Datasheet

Human MUC1, His Tag (orb1805342)

Description	Human MUC1, His Tag is expressed from HEK293. It contains AA 1135 (Accession No. P15941).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	MUC1
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Met 24 - Leu 1158
Isotype	His
Purity	> 85%
MW	123.2 kDa
Uniprot ID	P15941
Expiration Date	6 months from date of receipt.