

Product Datasheet

Human OLR1 / LOX1, His Tag (orb1805345)

Description	Human OLR1 / LOX1, His Tag is expressed from HEK293. It contains AA 216 (Accession No. P78380-1).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	OLR1 / LOX1
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Met 58 - Gln 273
Isotype	His
Purity	> 85%
MW	25.82 kDa
Uniprot ID	P78380
Expiration Date	6 months from date of receipt.