

Product Datasheet

Human Integrin b6 / ITGB6, Fc Tag (orb1805346)

Description	Human Integrin b6 / ITGB6, Fc Tag is expressed from HEK293. It contains AA 688 (Accession No. P18564...
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Integrin b6 / ITGB6
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Gly 22 - Pro 709
Isotype	Fc
Purity	> 85%
MW	201.16 kDa
Uniprot ID	P18564
Expiration Date	6 months from date of receipt.