

Product Datasheet

Human GPNMB, His Tag (orb1805348)

Description	Human GPNMB, His Tag is expressed from HEK293. It contains AA 476 (Accession No. Q14956).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	GPNMB
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Lys 23 - Asn 498
Isotype	His
Purity	> 85%
MW	54.4 kDa
Uniprot ID	Q14956
Expiration Date	6 months from date of receipt.