

Product Datasheet

Human Glycophorin A / CD235a, Fc Tag (orb1805350)

Description	Human Glycophorin A / CD235a, Fc Tag is expressed from HEK293. It contains AA 72 (Accession No. P027...
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Glycophorin A / CD235a
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	Fc
Note	For research use only
Protein Sequence	Ser 20 - Glu 91
Isotype	Fc
Purity	> 85%
MW	68.02 kDa
Uniprot ID	P02724
Expiration Date	6 months from date of receipt.