

## Product Datasheet

### Human Endoglin / CD105, His Tag (orb1805351)

<b>Description</b>	Human Endoglin / CD105, His Tag is expressed from HEK293. It contains AA 561 (Accession No. P17813?)...
<b>Species/Host</b>	Human
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 1 EU/mg
<b>Tested Applications</b>	ELISA, Kinetics
<b>Target</b>	Endoglin / CD105
<b>Preservatives</b>	1×PBS, pH 6.0
<b>Storage</b>	-80°C
<b>Tag</b>	His
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Glu 26 - Gly 586
<b>Isotype</b>	His
<b>Purity</b>	> 85%
<b>MW</b>	61.84 kDa
<b>Uniprot ID</b>	<b>P17813</b>
<b>Expiration Date</b>	6 months from date of receipt.