

Product Datasheet

Human Embigin, His Tag (orb1805353)

Description	Human Embigin, His Tag is expressed from HEK293. It contains AA 228 (Accession No. Q6PCB8).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	Embigin
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Asp 33 - Pro 260
Isotype	His
Purity	> 85%
MW	26.89 kDa
Uniprot ID	Q6PCB8
Expiration Date	6 months from date of receipt.