

Product Datasheet

Human TROP2, His Tag (orb1805363)

Description	Human TROP2, His Tag is expressed from HEK293. It contains AA 248 (Accession No. P09758).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	TROP2
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	His 27 - Thr 274
Isotype	His
Purity	> 85%
MW	28.98 kDa
Uniprot ID	P09758
Expiration Date	6 months from date of receipt.