

Product Datasheet

Human PVRIG, His Tag (orb1805366)

Description	Human PVRIG, His Tag is expressed from HEK293. It contains AA 131 (Accession No. Q6DKI7).
Species/Host	Human
Reactivity	Human
Endotoxins	< 1 EU/mg
Tested Applications	ELISA, Kinetics
Target	PVRIG
Preservatives	1×PBS, pH 6.0
Storage	-80°C
Tag	His
Note	For research use only
Protein Sequence	Thr 41 - Asp 171
Isotype	His
Purity	> 85%
MW	14.76 kDa
Uniprot ID	Q6DKI7
Expiration Date	6 months from date of receipt.